

Abicin 30% Resin Lacquer Questions and Answers

Question	Answer
What is Abicin 30% Resin Lacquer?	Abicin 30% Resin Lacquer is a modern CE Marked (EU Medical Device class IIb) lacquer for treatment of fungal nail infections, fungal skin infections (Athlete`s foot) and paronychia.
What are the indications for Abicin 30% Resin Lacquer?	Abicin 30% Resin Lacquer is indicated for home use in the treatment of of fungal nail infections, fungal skin infections (Athlete`s foot) and paronychia.
What are the main properties of Abicin?	The antifungal and antibacterial effectiveness of Abicin 30% Resin Lacquer has been studied and confirmed by <i>in vitro</i> and <i>in vivo</i> research methods required by modern medical science. Abicin 30% Resin Lacquer is effective against all the usual dermatophyte i.e. filamentous fungal infections of the skin and nails of man. The antifungal efficacy covers dermatophytes of the Trichophyton, Microsporum and Epidermophyton -geni. Abicin 30% Resin Lacquer is also very effective against many bacteria causing paronychia, such as Staphylococci and Streptococci.
How long should Abicin be used?	Abicin 30% Resin Lacquer is used for 3-12 months depending on the position and the severity of the fungal nail infection. The treatment of bacterial paronychia lasts 1-3 weeks.
Does Abicin 30% Resin Lacquer have an effect on inflammation?	Yes, Abicin 30% Resin Lacquer has strong anti-inflammatory properties.
Can diabetics use Abicin?	Yes, diabetic patients can use Abicin.
Is the product “sterile”?	Abicin is not labelled as being “sterile” in that it has not undergone any external sterilisation process. However the product is self-preserving due to its antimicrobial activity and is suitable for reuse over a 12 month period.
Is there bacterial resistance to Abicin?	No, there have not been any reports of bacterial resistance.
How is the resin cultivated?	When a tree is damaged it produces a sticky resin to plug the wound and protect the tree from infection. Over a period of 3 or more years this resin seals the tree wound. At this time the resin is harvested by hand.
How is the resin cleaned?	To purify the resin it must be separated from dirt and bark chippings and other unwanted components. This is done by a patented method at low temperature which prevents the production of potentially hazardous tars.
Can Abicin be stored in the refrigerator?	Yes. Abicin 30% Resin Lacquer can be stored in the refrigerator. However, it is not mandatory. It should not be stored in the freezer or in temperatures over 40°C.